

SIGMA Spectrum Infusion System and Spectrum IQ

Compatible Cleaners

- 10% Solution of Clorox® Bleach and 90% Water
- 3M™ HB Quat Disinfectant Cleaner Concentrate
- A-456 II Disinfectant Cleaner
- Accel Prevention Wipes
- Accel TB Wipes
- Clorox Commercial Solutions® Clorox® Germicidal Wipes
- Clorox Healthcare Hydrogen Peroxide Wipe
- Dispatch® Hospital Cleaner Disinfectant Towels with Bleach
- EnVerros Sanimaster 4®
- Oxivir TB Wipes
- Oxycide
- PDI® Sani-Cloth® Bleach Wipe
- Super Sani-Cloth® Germicidal Disposable Wipe
- Virex II 256
- Virox 5 RTU

Note:

- The listing is strictly to identify materials compatibility with the SIGMA Spectrum Infusion Pump and Spectrum IQ.
- Consult the respective Operator's Manual for instruction and cautions for cleaning.

This document is subject to change, please visit www.spectrumiq.com or contact Technical Support for information:

- **Phone: (800) 356-3454 and follow prompts**
- **Email: MedinaTechSupport@baxter.com**

Rx Only. For the safe and proper use of this device, refer to appropriate operator's manuals